



Specialty Lenses

Specialty Lenses

Available with the nanoCLEAR System

nanoCHROME and nanoFLASH lenses include a flash mirror front surface that gives the lens a mirror-like appearance

Available in six eye catching reflected colors

- Firestorm Red
- Desert Orange
- Electric Yellow
- Forest Green
- Ocean Blue
- Midnight Violet

nanoCHROME™

- Flash mirror coatings combined with your tint colors or photochromics provide nearly limitless fashion options

nanoFLASH™

- Untinted

nanoXT™

- UV-curable, nanoparticle hard coat
- Double-sided process
- Tintable

nanoGUARD™

- Super scratch resistant, dual layer treatment
- Performs like a factory dip hardcoat

Power Range

- Sphere +4.00D to -6.00D (0.25D steps)
- Cylinder -0.25D to -2.50D (0.25D steps)
- Add Power +1.00D to +3.00D (0.25D steps)

